



BostonSight Wet Shipping Notification

To aid with conditioning, each lens is supplied (non-sterile) in a plastic lens storage case filled with Boston SIMPLUS® Multi-Action Solution.

Boston SIMPLUS® Multi-Action Solution is a cleaning, disinfecting and conditioning product that has an added protein remover.

The solution needs to be replaced with fresh, sterile, and unexpired solution within 30 days from the manufacture date indicated on the device label. Refer to the Package Insert for additional instructions on lens care and handling.

Additional information on Boston SIMPLUS® Multi-Action Solution is available at www.bausch.com/ecp/our-products/contact-lens-care/gas-permeable-lens-care/boston-simplus-solution including the package insert which contains the solution composition as well as important safety and use information.

L-007 Rev 1